

## Instructions for Use

Version: 1.0.1  
Revision date: 13-Sep-24

### TUNEL In Situ Cell Apoptosis Detection Kit

Catalog No.	Dye
abx092371	FITC
abx092372	AF488
abx092373	AF594
abx092374	AF647
abx092375	AF555

**Size:** 50 tests / 100 tests

**Storage:** Store all components in the dark at -20 °C. Avoid repeated freeze/thaw cycles.

**Application:** For detecting cell apoptosis in tissue (paraffin-embedded or frozen) and cells (cell slides and cell smears).

#### Introduction

During apoptosis, specific DNA endonucleases are activated, cleaving genomic DNA between nucleosomes. The DNA of apoptotic cells is cleaved into multimers of 180-200 bp fragments. Fluorescent-labeled nucleotides can be attached to the fragments via the exposed 3'-OH end of the cleaved DNA, catalyzed by terminal deoxynucleotidyl transferase, which can be analyzed using a fluorescence microscope.

#### Kit Components (50 tests)

1. TdT Equilibration Buffer 9 ml
2. TdT Enzyme: 250 µl
3. Proteinase K (100X): 50 µl
4. Labelling Solution (see table above): 5 × 100 µl
5. DNase I (20 U/µl): 13 µl
6. DNase I Buffer (10X): 700 µl
7. DAPI Reagent (25 µg/ml): 250 µl

#### Kit Components (100 tests)

1. TdT Equilibration Buffer 2 × 9 ml
2. TdT Enzyme: 2 × 250 µl
3. Proteinase K (100X): 100 µl
4. Labelling Solution (see table above): 10 × 100 µl
5. DNase I (20 U/µl): 25 µl
6. DNase I Buffer (10X): 1500 µl
7. DAPI Reagent (25 µg/ml): 500 µl

#### Materials Required But Not Provided

1. Fluorescence microscope
2. PBS
3. Double-distilled water
4. Anti-fluorescence quenching reagent
5. Centrifuge and centrifuge tubes
6. Pipettes and pipette tips
7. Filter paper
8. Polyformaldehyde
9. Triton-100
10. Xylene
11. Ethanol

## Instructions for Use

Version: 1.0.1  
Revision date: 13-Sep-24

### Protocol

#### A. Reagent Preparation

##### 1. Proteinase K Working Solution

Dilute the 100X Proteinase K 100-fold with PBS to prepare the 1X Proteinase K Working Solution (e.g. add 1  $\mu$ l of 100X Proteinase K to 99  $\mu$ l of PBS to prepare 100  $\mu$ l of 1X Proteinase K Working Solution). Mix well. Prepare only prior to immediate use.

##### 2. DNase I Buffer Solution

Dilute the 10X DNase I Buffer 10-fold with double-distilled water to prepare the 1X DNase I Buffer Solution. Mix well. Prepare only prior to immediate use.

##### 3. DNase I Working Solution

Dilute the DNase I (20 U/ $\mu$ l = 20,000 U/ml) 100-fold with 1X DNase I Buffer Solution to prepare the DNase I Working Solution (0.2 U/ $\mu$ l = 200 U/ml). Mix well. Prepare only prior to immediate use. Do not vortex.

##### 4. DAPI Reagent Working Solution

Dilute the DAPI Reagent 25-fold with PBS to prepare the DAPI Reagent Working Solution (e.g. add 4  $\mu$ l of DAPI Reagent to 96  $\mu$ l of PBS to prepare 100  $\mu$ l of DAPI Reagent Working Solution). Mix well. Prepare only prior to immediate use.

##### 5. Fixative Buffer

Dissolve polyformaldehyde in PBS to a final concentration of 4%.

##### 6. Cell Permeabilization Buffer

Dissolve Triton-100 in PBS to a final concentration of 0.2%. Once prepared, the solution can be stored at 4 °C for up to 2 days.

##### 7. Labelling Working Solution

For each sample slide and positive control slide, prepare 50  $\mu$ l of Labelling Working Solution by mixing 35  $\mu$ l of TdT Equilibration Buffer, 10  $\mu$ l of Labelling Solution, and 5  $\mu$ l of TdT Enzyme. Prepare only prior to immediate use.

##### 8. Negative Control Labelling Working Solution

For each negative control slide, prepare 50  $\mu$ l of Negative Control Labelling Working Solution by mixing 40  $\mu$ l of TdT Equilibration Buffer and 10  $\mu$ l of Labelling Solution. This solution does not contain TdT Enzyme. Prepare only prior to immediate use.

#### B. Sample Fixation and Permeabilization

##### • Cell Samples

- 1. Cell slides:** Wash slides in PBS. Use filter paper to absorb excess liquid around the sample. Immerse slides in Fixative Buffer at room temperature for 15-20 minutes, or at 4 °C for 1-2 hours.

**Cell smears:** Collect cells and resuspend in PBS. Add an equal volume of Fixative Buffer to the volume of PBS previously added. Allow the mixture to stand at room temperature for 15-20 minutes, or at 4 °C for 1-2 hours. Centrifuge at 600  $\times$  g for 5 minutes. Discard the supernatant, then add PBS to resuspend the cells. On a slide, aliquot 25-50  $\mu$ l of cell suspension and allow to dry.

2. Wash slides with PBS 3 times, 5 minutes each time.
3. Immerse slides in Cell Permeabilization Buffer and incubate at 37 °C for 10 minutes.
4. Wash slides with PBS 3 times, 5 minutes each time.

## Instructions for Use

Version: 1.0.1  
Revision date: 13-Sep-24

### • Paraffin Sections

1. Deparaffinize and hydrate paraffin slides using conventional methods. Immerse slides in xylene 2 times, 10 minutes each time (may be increased to 20 minutes if the ambient temperature is lower than 20 °C). Hydrate paraffin sections with decreasing concentrations of ethanol (100%, 95%, 90%, 80%, 75%), 3 minutes each.
2. Wash slides with PBS 3 times, 5 minutes each time.
3. Use filter paper to absorb excess liquid around the sample. Add 100 µl of 1X Proteinase K Working Solution to each sample. Incubate at 37 °C for 20 minutes. The incubation time can be adjusted depending on the sample.
4. Wash slides with PBS 3 times, 5 minutes each time.

### • Frozen Sections

1. Equilibrate frozen sections to room temperature, then immerse frozen slides in Fixative Buffer at room temperature for 15-20 minutes.
2. Wash slides with PBS 2 times, 5 minutes each time.
3. Add 100 µl of 1X Proteinase K Working Solution to each sample. Incubate at 37 °C for 10-20 minutes. The incubation time can be adjusted depending on the sample.
4. Wash slides with PBS 3 times, 5 minutes each time.

### C. Positive and Negative Control Preparation

It is recommended to select representative samples as positive and negative controls. Positive controls are treated with DNase I to produce exposed 3'-OH ends; and are used to verify the effectiveness of the experimental process and reagents. Negative controls are treated with a labelling solution that does not include TdT Enzyme; and are used to determine non-specific staining of samples and sample autofluorescence.

1. Add 100 µl of 1X DNase I Buffer to each positive control slide and negative control slide. Allow the slides to stand at room temperature for 5 minutes.
2. **Positive control slides:** Use filter paper to absorb excess liquid around the sample. Add 100 µl of DNase I Working Solution (200 U/ml) to each positive control slide, then incubate at 37 °C for 10-30 minutes.

**Negative control slides:** Incubate at 37 °C for 10-30 minutes.

3. Wash slides with PBS 3 times, 5 minutes each time.

### D. Labelling

#### Notes:

- Bring the TdT Equilibration Buffer to room temperature and mix well before use. Crystals may be observed when the TdT Equilibration Buffer is taken out of the freezer; they will dissolve at room temperature. Ensure all crystals are dissolved prior to use.
- Store the TdT Enzyme at -20 °C immediately after the labelling working solutions are prepared. Do not vortex the TdT Enzyme.
- Before using the labelling working solutions, keep on ice and mix well.

#### Labelling Protocol

1. Add 100 µl of TdT Equilibration Buffer to each sample. Incubate at 37 °C for 10-30 minutes.
2. Use filter paper to absorb excess liquid around the sample. Do not allow the samples to dry out completely.

## Instructions for Use

Version: 1.0.1  
Revision date: 13-Sep-24

3. **Sample slides and positive control slides:** Add 50 µl of Labelling Working Solution to each sample slide and positive control slide.

**Negative control slides:** Add 50 µl of Negative Control Labelling Working Solution to each negative control slide.

4. Incubate at 37 °C for 1 hour in the dark in a humidified chamber. If the signal intensity is low, the incubation time can be extended up to 4 hours.

5. Wash slides with PBS 3 times, 5 minutes each time.

6. Use filter paper to absorb excess liquid around the sample. Add DAPI Reagent Working Solution and allow to stand in the dark for 5 minutes.

7. Use filter paper to absorb excess liquid around the sample. Add Anti-fluorescence quenching reagent to seal the slides.

### E. Analysis

Samples can be analyzed using a fluorescence microscope with an appropriate filter. It is recommended to analyze samples as soon as possible or store samples at 4 °C in the dark.

Dye	Ex/Em (nm)	Filter Set
FITC	490/520	FITC Filter Set
AF488	495/519	FITC Filter Set
AF594	590/617	TRITC Filter Set
AF647	650/665	Cy5 Filter Set
AF555	555/565	TRITC Filter Set
DAPI	350/470	DAPI Filter Set

### F. Technical Support

For troubleshooting and technical assistance, please contact us at [support@abbexa.com](mailto:support@abbexa.com).