

Beagle dog Whole Blood (Sterile Sodium Citrate)

Catalogue No.: abx800102

This product is currently in development. The lead time for this product may be several months. Please contact us at info@abbexa.com for an updated lead time before purchasing this product.

Whole blood fraction from Beagle dog in Sodium Citrate. This product is Sterile. Can be used for titration of complement, adsorption procedures, testing for agglutinins/HA assays, and for the preparation of stroma as particulate reagents.

Target:	Beagle dog Whole Blood (Sterile Sodium Citrate)
Origin:	Beagle dog
Recommended dilutions:	Optimal dilutions/concentrations should be determined by the end user.
Form:	Liquid
Purification:	Beagle dog whole blood was harvested without any purification steps. This product was aseptically drawn and processed using sterile equipment and sterile techniques.
Storage:	Store at 4°C.
Validity:	The exact shelf may depend on the application. Whole blood is a perishable product and can be preserved for up to 6 weeks.
Buffer:	Sodium Citrate.
Note:	THIS PRODUCT IS FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC, THERAPEUTIC OR COSMETIC PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION.
Directions for use:	Use aseptic techniques to maintain sterility when opening the product.