## **Datasheet**

Version: 2.0.0 Revision date: 06 Oct 2024



## Low Sample Volume Lipopolysaccharides (LPS) ELISA Kit

Catalogue No.:abx052841

Lipopolysaccharides (LPS) ELISA Kit is an ELISA Kit for the in vitro quantitative measurement of Lipopolysaccharide (LPS) concentrations in serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids. This kit is specifically designed for the detection of low volume samples.

Target: Lipopolysaccharides (LPS)

Reactivity: General

Tested Applications: ELISA

Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.

Storage: Shipped at 4 °C. Upon receipt, store the kit according to the storage instruction in the kit's manual.

Validity: The validity for this kit is at least 6 months. Up to 12 months validity can be provided on request.

**Stability:** The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5% within

the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested

that the whole assay is performed by the same user throughout.

**Test Range:** 12.35 ng/ml - 1000 ng/ml

Sensitivity: < 4.53 ng/ml

Standard Form: Lyophilized

**Detection Method:** Colorimetric

Assay Type: Competitive

Assay Data: Quantitative

Sample Type: Serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological

fluids.

## **Datasheet**

Version: 2.0.0 Revision date: 06 Oct 2024



Note:

This product is for research use only.

The range and sensitivity is subject to change. Please contact us for the latest product information. For accurate results, sample concentrations must be diluted to mid-range of the kit. If you require a specific range, please contact us in advance or write your request in your order comments. Please note that our ELISA and CLIA kits are optimised for detection of native samples, rather than recombinant proteins. We are unable to guarantee detection of recombinant proteins, as they may have different sequences or tertiary structures to the native protein.



Website: www.abbexa.com · Email: info@abbexa.com