

## Human Metalloproteinase Inhibitor 2 (TIMP2) Protein (Active)

Catalogue No.: abx655837

Human Tissue Inhibitors Of Metalloproteinase 2 (TIMP2) Protein (Active) is an Active Recombinant Human protein.

<b>Target:</b>	Metalloproteinase Inhibitor 2 (TIMP2)
<b>Origin:</b>	Human
<b>Expression:</b>	Recombinant
<b>Tested Applications:</b>	WB, SDS-PAGE
<b>Conjugation:</b>	Unconjugated
<b>Form:</b>	Lyophilized
<b>Purity:</b>	> 90%
<b>Reconstitution:</b>	To keep the original salt concentration, we recommend reconstituting to the original concentration prior to lyophilization (see Concentration) in ddH <sub>2</sub> O. If a lower concentration is required, dilute in 20 mM Tris, 150 mM NaCl, pH 8.0. If a higher concentration is required, the product can be reconstituted directly in 20 mM Tris, 150 mM NaCl, pH 8.0, though please note that this will change the overall salt concentration. The stock concentration should be between 0.1-1.0 mg/ml. Do not vortex.
<b>Storage:</b>	Store at 2-8 °C for up to one month. Store at -80 °C for up to one year. Avoid repeated freeze/thaw cycles.
<b>Sequence Fragment:</b>	Please enquire.
<b>Tag:</b>	N-terminal His tag
<b>Buffer:</b>	Prior to lyophilization: 20 mM Tris, 150 mM NaCl, pH 8.0, containing 1 mM EDTA, 1 mM DTT, 0.01% Sarcosyl, 5% Trehalose and Proclin-300.
<b>Activity:</b>	Active
<b>Concentration:</b>	Prior to lyophilization: 200 µg/ml
<b>Note:</b>	This product is for research use only. Not for human consumption, cosmetic, therapeutic or diagnostic use.