

# Datasheet

Version: 5.0.0  
Revision date: 03 Apr 2025



## Protein Loading Buffer (6X)

Catalogue No.: abx098124

6X Protein Loading Buffer for preparing protein samples for SDS-PAGE. Prior to loading, add 6X Protein Loading Buffer to protein sample to a working concentration of 1X.

**Target:** Protein Loading Buffer

**Storage:** Store at -20°C. Avoid repeated freeze/thaw cycles.

**Validity:** Up to 2 years.

**Buffer:** The exact formulation is proprietary.

**Concentration:** 6X

**Note:** THIS PRODUCT IS FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC, THERAPEUTIC OR COSMETIC PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION.

For Reference Only